




- Silicone adhesive backing



- Foam Layer

Conlett Neo Guard Plus is a further development of our Neo-Guard series of products. It is designed to protect sensitive neonatal nasal and facial skin during respiratory therapy.

With a layer of foam and latest technology of silicone dressing, it will effectively protect the neonatal septum from trauma and irritation caused by pressure and motion of prongs.



- Star Cut Hole

Features of the Neo Guard Plus:

- Silicone is super skin friendly
- There is an extra cushion layer of foam to protect the nasal septum
- Repositioning of prongs is possible
- A cut-through star around the holes to generate retention and create a tight seal

SIZES AVAILABLE

When ordering please quote the following:

DESCRIPTION	ORDER#	QTY/BOX
CPAP – Neo Guard Plus w/foam (Small to Medium) <2000g	NGP-S-M	50
CPAP – Neo Guard Plus w/foam (Medium to Large) >2000g	NGP-M-L	50

(Minimum order quantity of 1 box)

Conlett

U25/17 Cairns St
Loganholme
QLD 4129
Australia

Tel: (+61) 7 3416 8902
Email: contact@conlett.com