

WaltersMedical

Specialist Critical Care



Neonatal Products

Consumables & Devices

www.waltersmedical.co.uk

Atom Products

To save a tiny baby's life

Seeing really is *believing*...

**Proud distributor of Atom
Incubators & Associated Products***

The Atom range includes:

- Incubators
- Infant Warmers
- Transport Incubators
- Phototherapy
- Bolt-On Accessories

* In England, Scotland & Wales only



Dual Incu i

CAT NO. AT61601

This model can perform in a wide range of hospital settings through the incubator mode and infant warmer mode.



Please visit our website for more details on the Atom range.

www.waltersmedical.co.uk

A thriving partnership



Available via
**NHS Supply
Chain**

**BOOK
YOUR
DEMO**



Infra Warmer i

CAT NO. AT61634

High-end Infant warmer with various features enabling its use in either acute or chronic environments.



Jaundice & Phototherapy

Neonatal Transcutaneous Jaundice Meter

CAT NO. BM-100A

Designed for departmental and community use and securely protected by a metal carry case, our new device provides a reliable solution for measuring bilirubin levels in newborns.

Annually Calibrated & Serviced at our UK Service Centre



Rechargeable
battery capable
of performing
2,000
measurements

**REQUEST
A DEMO**



bili-hut™

CAT NO. 0130A

**Revolutionising Neonatal
Phototherapy at Home**

The bili-hut™ brings hospital-grade phototherapy care to families, reducing inpatient stays and readmissions.

Perfect for neonatal outreach services as well as Transitional Care & Postnatal Ward use.

Neo-Photoshades

Provide a protective and secure eye shield for the infant when they require phototherapy treatment; an ergonomic shape, soft cushioned structure and extremely gentle silicone adhesive ensures the Neo-Photoshades block any harmful light without damaging sensitive or fragile neonatal skin.

CAT NO.	SIZE	WEIGHT	BOX QTY
NSPS000	Petite	≤900g	50
NSPS001	Small	<1000g	50
NSPS002	Medium	>1000g	50



Available
via the
NHS Supply
Chain

Neo Positioning

IV Armboard with Straps

Manufactured from soft polyurethane foam with an aluminium core strip and soft straps with hook & loop fastening for gentle, comfortable and secure immobilization of the neonate's limb.

- Simple application, repositioning and removal
- Malleable soft foam cover for added comfort and protection
- Aluminium core ensures profile maintained
- Gentle, soft fastening straps designed for ease of use
- Available in 4 sizes



CAT NO.	NAME	NHSSC NPC	BOX QTY
NBS001	Armboard with straps 70 x 29mm (3")	FVH85250	50
NBS002	Armboard with straps 110 x 41mm (4")	FVH85249	50
NBS003	Armboard with straps 139 x 50mm (5")	FVH85248	50
NBS004	Armboard with straps 175 x 50mm (7")	FVH85247	50



Neo-Hug

An infant positioning aid providing optimum support boundaries whilst offering a secure, warm, environment to the cocooned neonate. Available in two sizes with disposable covers.

CAT NO.	SIZE	LENGTH	HEIGHT	BOX QTY
NH001	Petite	90cm	6.5cm	5
NH002	Small	90cm	9cm	5

Neo-Teneo

A non-toxic, gell-filled, positional cushion to help alleviate pressure caused by prolonged immobilisation where frequent turning or repositioning of the neonate is not possible. Use of Neo-Teneo can help to reduce Plagiocephaly and absorption of vibrations whilst protecting fragile preterm skin.

- Soft and gentle on neonate's skin
- Reusable and easy to clean

CAT NO.	SHAPE	DIMENSIONS	BOX QTY
NT-001	Rectangular	13 x 19 x 2cms	15



Disposable covers are also available; please contact us for details


Neo Protection

Neo-Guard and Neo-Guard Plus

Protects sensitive neonatal nasal and facial skin during respiratory therapies, effectively guarding the neonatal septum from trauma and irritation caused by continual pressure and motion of prongs and masks.

- Super skin friendly and gentle adhesive for sensitive, delicate, neonate skin
- Repositionable
- Easy observation of the condition of the nasal nares
- Multiple sizes & designs for the best combination of therapies and baby weights
- **Neo-Guard Plus** integrates an extra layer of foam to further protect the nasal septum


CAT NO.	TYPE	WEIGHT	BOX QTY
NG000	CPAP Neo-Guard	< 700g	50
NG001	CPAP Neo-Guard	700g - 1250g	50
NG002	CPAP Neo-Guard	1250g - 2000g	50
NG003	CPAP Neo-Guard	2000g - 3000g	50
NG004	CPAP Neo-Guard	> 3000g	50
NGP-S-M	CPAP Neo-Guard Plus w/foam	< 2000g	50
NGP-M-L	CPAP Neo-Guard Plus w/foam	> 2000g	50



Neo-Guard Plus
Please contact us for further details or to discuss your specific requirements

Neo-Flexi Thermal Wrap

Helps to maintain the temperature of newborn, preterm and low birth-weight babies, by preventing heat loss and hypothermia through evaporation, convection and conduction. Suitable for use during resuscitation, stabilization and for transfer/transport of newborns to specialised neonatal intensive care units. One size only.



CAT NO.	SIZE	BOX QTY
FTW-001	60cm x 60cm	50

Sonic Ear Muffs

Protecting the sensitive cochlear of premature babies, our muffs are ergonomically designed to fit comfortably around the baby's ear, shielding background noise disturbance to facilitate a quiet and calm developmental environment.

Useful during MRI and transport of infants



30 PAIRS PER BOX

CAT NO. SM-001

Neo Treatment Aids

Disposable SpO2 Probes

We offer a complete range of disposable and reusable wrap sensors, including Oximax compatible sensors.

- Choose from gentle adhesive or non-adhesive extra soft foam

Please contact us for further details or to discuss your specific requirements.

NO.	CAT NO.	NHSSC NPC	WEIGHT
1	SSD-001XC-W09AN	FBC85035	Neo/Adult up to 40kg
2	SSD-001XC-W09I	FBC85034	Infant 3-20kg
3	SSD-001XC-W09P	FBC85037	Paediatric 10-50kg
4	SSD-001XC-S09PI	FBC85032	Infant/Paed 3-40kg
5	SSD-001XC-S09AN	FBC85033	Adult /Neo >40kg



Available
via the
NHS Supply
Chain

Neonatal/Paediatric Foam Wraps & Moveable ID Bracelet

Manufactured from an extremely soft foam, and with a choice of three variants, our wraps snugly secure all types of SpO2 sensor in place eliminating the use of harsh adhesive tape to preserve the delicate skin of the baby/patient.

Just like our soft foam wraps, our Moveable ID Bracelet is gentle on the neonate's delicate skin but includes a writing panel as an added benefit.

CAT NO.	NHSSC NPC	SIZE	WIDTH	LENGTH	BOX QTY
WMW6664H	FBC1591	Large with holes	3cm	12.5cm	12
WMW6554	FBC1603	Large	3cm	12.5cm	12
WMW6444S	FBC1602	Small	3cm	9cm	15
WMW4648	FBC1590	Large with ID Panel	3cm	12.5cm	12



NEONATAL BEST SELLER WMW6554

WE ALSO OFFER A FULL RANGE OF NIBP HOSES TO ACCOMPANY THESE CUFFS

Available
via the
NHS Supply
Chain

Neonatal NIBP cuffs

Our cuffs are manufactured to a very high standard to ensure durability and comfort for the baby without compromising on price.

- Made from extremely soft PU material
- Available in 5 neonate sizes
- Option of single or double tube cuff
- Colour-coded cuffs & boxes for fast size identification



Neo-Vista & Neo-Vista Plus Transilluminators

For easy identification of veins in neonates and babies our Neo-Vista portable illumination devices illuminate a map of the vasculature on the surface of the skin aiding all needle procedures.

- Enables clinicians to verify patency and avoid valves or bifurcations
- Increases success rate with first attempt of venepuncture



Sil-Fixx Tape

A transparent, soft and extremely pliable fixation tape especially suitable for very fragile and sensitive neonatal skin.

- Excellent adhesive properties leaving no residue
- Can be re-applied without affecting the adhesive quality
- Does not stick to hair / Waterproof
- All types of fixation e.g. IV, Cannula, Tubes, Electrodes etc



Neo-Shield Skin Temperature Probe Covers

The latest technology against heat radiation and convection to ensure accurate temperature monitoring.

Extremely gentle Silicone based adhesive to maintain the neonate's skin integrity; will not cause excoriation when removed.

CAT NO.	SIZE	BOX QTY
REF-001	2.5 x 2.5cm	50
REF-002	3.5 x 3.5cm	50



Disposable Temperature Probes

CAT NO.	FOR USE WITH	BOX QTY
TSD-DGB-GS30 W0026E	Draeger Neonatal Devices GE Neonatal Devices	25 24

CAT NO.	SIZE	PACK SIZE
SF1915	1.9cm x 1.5m roll	6 rolls
SF1950	1.9cm x 5m roll	6 rolls
SF2515	2.5cm x 1.5m roll	6 rolls
SF2550	2.5cm x 5m roll	6 rolls
SF5015	5cm x 1.5m roll	4 rolls
SF5050	5cm x 5m roll	4 rolls



We're here to help...

REQUEST A
SAMPLE BOX
TODAY



Caring for newborns is at the heart of everything we do.

Our neonatal consumables and devices are designed with comfort, safety, and simplicity in mind. If we can support your team's care and procurement needs, please call or email to arrange a meeting, or request a complimentary sample box to experience the quality for yourself.

Do you have a Walter Teddy Bear?

Share his journey on socials with us #WalterTheBear

Don't forget... Discover our Labels

With a catalogue including many popular neonatal-specific syringe and line label titles, as well as a bespoke label service catering for maternity, CTG and milk kitchen labels be sure to contact us for your label requirements

**ALSO
AVAILABLE**

Wall mounted &
free standing
Perspex
containers



500
LABELS
PER BOX



REQUEST A
FREE LABEL
SAMPLE
PACK

Please contact us to discuss your requirements or follow us on socials:

t: 0845 6800 377 e: sales@waltersmedical.co.uk

www.waltersmedical.co.uk

