WaltersMedical

Specialist Critical Care









Vein Finder

CAT NO. WMQV600



REQUEST A DEMO TODAY

www.waltersmedical.co.uk

WMQV600 Vein Finder

Introducing the WMQV600 Vein Finder.

Advanced vision for fast, confident, accurate access.

The WMQV600 is a portable vein finding device using infrared light to detect veins beneath the skin, projecting the position of the veins onto the skin surface directly above.

Designed to assist medical staff in precise vein location in complex situations - ***first time!***

Benefits include reducing patient pain and anxiety, saving time and resources and minimising complications.

The ergonomic design ensures the device is comfortable to hold, more flexible to use and more convenient to operate.

Features

The WMQV600 Vein Finder has various user-defined functions that can adapt to different ages, skin tone, weights and various operating environments.

- **Colour –** there are a choice of 12 colours – subject to user preference and adjustable to perform effectively across all skin tones
- **Size –** there are three image sizes suitable for infant, child and adult
- **Inverse –** this setting reverses the graphics to an alternative view
- **Brightness –** adjustable to your environment
- **Mode –** includes options to increase clarity of smaller, deeper veins, remove possible hair interference and locate the best puncture position
- **Photograph –** this enables extraction of a JPEG image from the device using the USB-C cable. 6,000 image capacity.




Instructional video available














Specification

Distance	15-25cm
Infrared Wave	760-940nm
Noise	≤30dB
Weight	450g
Product Size	24cm(L) 6cm(W) 6.5cm(H)
Battery	4900mAh rechargeable lithium battery
Power	DC 5V 4.5A
Adaptor	100V-240V, 50Hz-60Hz AC, 5V 2A
LCD Screen	2.4 inch
Image	3D pseudo-colour and AI algorithm
Accuracy of Blood Vessel Position	±0.3mm
Infrared Detection Depth	0-10mm

For hands free use the WMQV600 has the option of a mobile or table top stand.

Mobile Stand

The mobile stand is the most popular choice allowing the device to be securely moved between patients.

CAT NO. WMQV600-MS

Table Stand

The table stand is perfect for use in one specific location for example an Outpatients Clinic Room.

CAT NO. WMQV600-TS



Widely used in the following departments:

- **Elderly Care**
- **Emergency**
- **Laboratory**
- **Oncology**
- **Outpatients**
- **Phlebotomy**
- **Plastic Surgery**
- **Radiology**
- **Surgery**
- **Vascular**
- **Paediatrics** *and more*

WMQV600 Standard Configuration Includes:

• Power Adaptor • USB-C Cable • Tourniquet • Skin Marker • User Manual • Carry Case

We're here to *help*...

Walters Medical was born from our extensive knowledge and experience in medical electronics and consumables within hospital environments. We understand the importance of supporting patient care, clinical efficiency, and daily hospital operations, and are committed to supplying quality products without compromising on price. Please contact us if we can assist you.

WMQV600 supplied with a carry case



After Care Service

We offer full repairs and servicing for the WMQV600 at our UK Service Centre. Contact us for details.


You may also be interested in...

Neo Vista Transilluminators

Portable vein illumination devices for neonatal and paediatric use.



Neo-Vista TRANS1



Neo-Vista Plus TRANS2

Specialist Supplier:

ECG Cables • SpO2 & Temperature Probes • NIBP & IBP Products • Medical Labels
Neonatal & Protection Products • Patient Assessment Devices • Atom Incubators & Accessories • ***Repair Service & Parts***

Trusted expertise to find the product you need, fast.

www.waltersmedical.co.uk

WaltersMedical
Specialist Critical Care

t: 0845 6800 377
e: sales@waltersmedical.co.uk

Walters Medical Ltd 30 Campus 5, Third Avenue, Letchworth Garden City, Herts, SG6 2JF, UK

Keep up to date

Follow us on Socials:

17317v2.1026